


Connect ScanCast™ or Barcode Scanner to Omnicell CPM

1. **Login** to Omnicell CPM
2. **Click** the **“Gear Icon”** (for Settings) on the left-hand side of the Omnicell CPM top level screen
3. **Click** the **“Scanner Utility”** tab on the top of the Settings screen
4. **Click** **“Acquire Scanner”** Button (Step 2 on screen)


Using Omnicell Scanner

5. **Scan** the pairing barcode



Use ConsortiEX Scanner

5. **Connect** ScanCast™ USB Cable to the Omnicell CPM Computer
6. **Scan** ScanCast™ CPM Reset Barcode



7. **Scan** the Pairing Barcode



“Scanner is ready to Use” popup is displayed on Omnicell screen to show pairing was successful.

Troubleshooting:

If the barcode fails to pair (no “Scanner is ready to use” popup) try using the Toggle LE (Line Ending) barcode and follow the above sequence again.

If the ScanCast™ acquires in CPM (“Scanner is ready to Use”), however, 2D barcodes do not cause any action by CPM; try scanning the Toggle GS barcode.

ScanCast™ will remember the settings of these barcodes, even over power cycles.

Toggle GS (Group Separator)



Toggle LE (Line Ending)

