



ScanCast™ Action Barcodes:


[BD] [BD Pharmogistics] [SCRXBDR] Medium (1.75x1.75)




Quantity




Confirm




Send to Queue



Reset



Enable ^[1]




Disable

[1] Default setting after power-on.


[2] This setting is non-volatile and will be retained over power-cycles.

ScanCast™ Action Barcodes:


[BD] [BD Pharmogistics] [SCRXBDR] Medium (1.75x1.75)



Line Ending Carriage Return ^[2]



Line Ending Line Feed ^[2]



Line Ending Carriage Return/
Line Feed ^[1,2]

[1] Default setting after power-on.

[2] This setting is non-volatile and will be retained over power-cycles.