







Medication Receiving Process: ScanCast™ MC020085 with Verify on Receipt™ Troubleshooting

Issue	Resolution
<p>Medication did not scan in relationship with the PO – the incorrect PO was active on the screen.</p>	<ol style="list-style-type: none"> 1. Scan “Undo Line Item”  2. Scan “Reset”  3. Choose the correct PO <ol style="list-style-type: none"> a. Continue with normal scanning process
<p>If a medication needs to be put somewhere that did not have a shelf location.</p>	<ol style="list-style-type: none"> 1. Scan the 1st item then manually enter where the medication is going ex. Another location in the pharmacy, clinic, etc.) 2. Scan “Skip Shelf” 
<p>Duplicate scan occurs - mistakenly scanned the same medication barcode twice. Upon ConsortiEX popup "Do you want to accept duplicate?"</p>	<ol style="list-style-type: none"> 1. Scan “Cancel Dialog”  2. Scan correct item to meet quantity expectations 3. Continue with normal scanning process

<p>In the event of mis-scanning or completely losing your place in the scanning process</p>	<p>1. Scan “Reset”</p> 
<p>2D barcode will not scan – unreadable or damaged</p>	<p>1. Contact Wholesaler/Supplier 2. Start the process to ship the medication back for replacement/credit</p>
<p>No output from either ScanCast™ - applications are not populating correctly. No flashing lights on the ScanCast™ unit (look for orange and green lights)</p>	<p>1. Unplug the USB power. 2. Wait 5 seconds. 3. Plug USB power back into port. 4. Check to see that lights are flashing. 5. Continue with normal scanning process.</p>
<p>2D barcode not in the database of a connected application (EHR, inventory system, automation, etc.)</p>	<p>1. Scan “Skip Next Barcode”</p>  <p>2. Scan 1D barcode:</p> <ol style="list-style-type: none"> If the 1D barcode is recognized, contact IT to have them enter the 2D information for the current product in the database If the 1D barcode is not recognized, set the product aside, contact IT and have them enter the entire new product into the database.

**After scanning the 2D barcode and receiving
the error “NDC not recognized”**

1. Scan “Undo Line Item”



2. Scan 1D barcode
 - a. If the 1D barcode fails, send photos of the following to Support@consortiox.com
 - i. Container label
 - ii. 2D barcode
 - iii. 1D barcode
 - iv. Medication package/box
 - v. In app screenshot